

500 µL, Model 1001.5 TLL SYR, Long Life Instrument Syringe

Specification Sheet (Part/REF # 8250-35)



Product Specifications

Volume	500 µL
Syringe Type	Longlife Gastight Half-Height
Syringe Series	1000
Termination	PTFE Luer Lock (TLL)
Plunger Style	X-Style
Fluid Path	Borosilicate glass, UHMWPE, PFA
Autoclavable	No
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.181 inches (4.61 mm)
Barrel Outer Diameter	0.355 inches (9.02 mm)
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/oem/syringes/8250-35>

Spec. Version B