

Validation / Calibration Kit for Incyte

Specification Sheet (Part/REF # 243740)



The Validation / Calibration Kit for Incyte contains all material needed for two calibrations and several calibration validation checks for Incyte Sensors, including the permittivity simulator (Ref. 243743), Solution A (Ref. 238988) and Solution B (Ref. 243742).

Product Specifications

Content Solution B, Standard conductivity solution 12,88 mS, Calibrator to check the internal electronics of the Incyte Pre-Amplifier

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/process-analytics/sensors/243740>