

pH Cleaning Solution Set

Specification Sheet (Part/REF # 238290)



A two-part cleaning solution for all glass membrane pH sensors. Also includes storage solution. Glass cleaning Solution B contains material of animal origin.

Product Specifications

Electrolyte	Cleaning Solution Set (2x500mL A+B + 1x500mL Storage Solution)
Package	2 x 500 mL A+B + 1 x 500 mL Storage Solution
Stability Months (after date of manufacturing)	24
Storage Conditions	5 - 35 °C / 41 - 95 °F

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/process-analytics/sensors/238290>

Spec. Version C