

DuraCal pH buffer 4.01/7.00/9.21

500 ml

Specification Sheet (Part/REF # 238922)



Three buffer set for calibration at common pH levels. Each liquid solution comes with certificate and serves as a Certified Reference Material.

Product Specifications

pH Accuracy	± 0.01 / ± 0.02 pH
Package	500 mL, each
Stability Months (after date of manufacturing)	± 0.01 pH: 24 / ± 0.02 pH: 60
Certificate	DAkkS
Storage Conditions	5 - 35 °C / 41 - 95 °F
pH Value	4.01/7.00/9.21

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/process-analytics/sensors/238922>

Spec. Version C