

Validation / Calibration Kit for Incyte LC

Specification Sheet (Part/REF # 243713)



The Validation / Calibration Kit for Incyte LC contains all material needed for two calibrations and several calibration validation checks for Incyte Sensors, including the permittivity simulator (Ref. 243743), Solution A (Ref. 238988) and Solution B LC (Ref. 243708).

Product Specifications

Content Solution B, Standard conductivity solution 12,88 mS, Calibrator to check the internal electronics of the Incyte Pre-Amplifier

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/process-analytics/sensors/243713>

Spec. Version B