

# Hamilton Electrolyte 3 M KCl, 500 mL

Specification Sheet (Part/REF # 238936)



Genuine Hamilton electrolyte to keep your sensor running smoothly. 500mL bottle containing Hamilton's 3M KCl electrolyte for pH and ORP Sensors

## Product Specifications

Electrolyte	3 M KCl
Package	500 mL
Stability Months (after date of manufacturing)	24
Storage Conditions	5 - 35 °C / 41 - 95 °F

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/process-analytics/sensors/238936>

Spec. Version C