

Calibration Station

Specification Sheet (Part/REF # 243575)



Calibration station used for validation / calibration of sensors in gaseous media. Equipped with two inlets for a wide array of test scenarios.

Product Specifications

For Sensors with	120 mm, 12 mm, PG13,5
Dimensions (W x D x H)	137 x 115 x 186 mm
Housing Material	Black PEEK with aluminium holder
Ø O.D. Fitting	8,7 mm
Ø I.D. Fitting	3,5 mm
Pressure Range bar g	0 to 12 bar

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/process-analytics/sensors/243575>

Spec. Version B