

Hamilton Electrolyte 3 M KCl viscous, 500 mL

Specification Sheet (Part/REF # 238939)



Replacement electrolyte for use in refillable lab sensors.

Product Specifications

Electrolyte	Viscous 3 M KCl-LR
Package	500 mL
Stability Months (after date of manufacturing)	24
Storage Conditions	5 - 35 °C / 41 - 95 °F

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/process-analytics/sensors/238939>

Spec. Version C