

Skylyte, 100 mL

Specification Sheet (Part/REF # 238037)



Ensure smooth operation of your Hamilton pH sensors with the Genuine Hamilton Skylyte Electrolyte, conveniently housed in a 100mL bottle containing Hamilton's Skylyte Electrolyte. Complete the maintenance on your sensor by using Hamilton cleaning and storage solutions to ensure high accuracy and long life.

Product Specifications

Electrolyte	Skylyte
Package	100 mL
Stability Months (after date of manufacturing)	24
Storage Conditions	5 - 35 °C / 41 - 95 °F

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/process-analytics/sensors/238037>

Spec. Version C