

Conductivity standard 12880 µS/cm, Basic Line

Specification Sheet (Part/REF # 238988)



Hamilton conductivity standard - 12880 µS/cm, 500mL, plastic bottle

Product Specifications

Accuracy at 25 °C	± 1 %
Value at 25 °C	12,880 µS/cm
Package	500 mL Calpack bottle
Stability Months (after date of manufacturing)	18
Certificate	DFM
Storage Conditions	5 - 35 °C / 41 - 95 °F

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/process-analytics/sensors/238988>

Spec. Version C