

# Conductivity standard 5 $\mu\text{S}/\text{cm}$

Specification Sheet (Part/REF # 238926)



Hamilton Conductivity standard - 5  $\mu\text{S}/\text{cm}$  300 mL

## Product Specifications

Accuracy at 25 °C	$\pm 1 \%$
Value at 25 °C	5 $\mu\text{S}/\text{cm}$
Package	300 mL glass bottle
Stability Months (after date of manufacturing)	36
Certificate	DFM
Storage Conditions	5 - 35 °C / 41 - 95 °F

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/process-analytics/sensors/238926>

Spec. Version C