

DuraCal pH buffer 1.09, 500 ml

Specification Sheet (Part/REF # 238271)



Patented buffer solutions provide pH stability of up to five years. High buffering capacity provides rapid and stable calibration.

Product Specifications

pH Accuracy	± 0.02 pH
Package	500 mL
Stability Months (after date of manufacturing)	60
Certificate	Hamilton
Storage Conditions	5 - 35 °C / 41 - 95 °F
pH Value	1.09

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/process-analytics/sensors/238271>

Spec. Version C