

# DuraCal pH buffer 6.00, 500 ml

Specification Sheet (Part/REF # 238276)



Patented buffer solutions provide pH stability of up to five years. High buffering capacity provides rapid and stable calibration.

## Product Specifications

<b>pH Accuracy</b>	± 0.02 pH
<b>Package</b>	500 mL
<b>Stability Months (after date of manufacturing)</b>	60
<b>Certificate</b>	Hamilton
<b>Storage Conditions</b>	5 - 35 °C / 41 - 95 °F
<b>pH Value</b>	6

Specifications are subject to change without notice

**Spec. Version C**

<https://www.hamiltoncompany.com/process-analytics/sensors/238276>