

Beverly (discontinued)

Specification Sheet (Part/REF # 817100)



The Beverly portable dissolved oxygen monitor provides a battery operated system for sample line measurements found in food and brewery applications.

Product Specifications

Measurement Principle	Oxygen dependent luminescence quenching, temperature compensation
Measuring Range	20 ppb to 25 ppm (DO)
Dimensions (W x D x H)	222 x 142 x 322 mm
Operating Temperature Range	0 to 80 °C
Pressure Range bar g	0 to 10 bar
Protection Rating	IP 67

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/process-analytics/sensors/817100>

Spec. Version B