

**10 µL, Model 701 N Shimadzu
Syringe, (metal flange), Cemented
Needle, 23s ga, 1.69 in, point style
AS**

Specification Sheet (Part/REF # 202643)



Product Specifications

Volume	10 µL
Syringe Type	Microliter
Syringe Series	700
Termination	Cemented Needle (N)
Needle	Cemented Needle
Gauge	23s gauge
Point Style	AS: Conical, Autosampler
Needle Length	1.69 in (42.9 mm)
Fluid Path	SST, borosilicate glass
Autoclavable	No
Gas Sterilizable	Yes
Maximum Temperature	50°C / 122°F
Minimum Temperature	68°F (20°C)
Needle Inner Diameter	0.0045 inches (0.11 mm)
Needle Outer Diameter	0.0238 inches (0.64 mm)
Std Needle Dead Volume	0.46 µL
Sales Quantity	1 syringe per pack

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/syringes/202643>

Spec. Version B