

Analytical Guard Column, PRP-1

Specification Sheet (Part/REF # 79286)



Protect your PRP-1 column.

Product Specifications

Separation Mode	Reversed-Phase
Stationary Phase	PRP-1
Format	Guard Column
Hardware Inner Diameter	2.1 mm
Hardware Length	33 mm
Particle Size	10 µm
USP Classification	USP L21
Packing Material Type	PSDVB
Particle Surface Area (m2/cm3)	400
Pore Volume (cm3/g)	0.772
Guard Column Size	Analytical
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
UOM	1 EACH

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/hplc-columns/79286>