

PRP-X100 Anion Exchange Analytical Guard Column

Specification Sheet (Part/REF # 79287)



Protect your PRP-X100 column.

Product Specifications

Separation Mode	Anion Exchange
Stationary Phase	PRP-X100
Format	Guard Column
Hardware Inner Diameter	2.1 mm
Hardware Length	33 mm
Particle Size	10 µm
Packing Material Type	PSDVB/Trimethylammonium
Particle Surface Area (m2/cm3)	405
Pore Volume (cm3/g)	0.64
Guard Column Size	Analytical
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
UOM	1 EACH

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/hplc-columns/79287>