

PRP-X300 Analytical Guard Column

Specification Sheet (Part/REF # 79289)



Organic acids and alcohols.

Product Specifications

Separation Mode	Ion Exclusion
Stationary Phase	PRP-X300
Format	Guard Column
Hardware Inner Diameter	2.1 mm
Hardware Length	33 mm
Particle Size	12-20 µm
Packing Material Type	PSDVB/Sulfonic Acid
Particle Surface Area (m ² /cm ³)	360
Pore Volume (cm ³ /g)	0.748
Guard Column Size	Analytical
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Compound Classification	Organic Acids/Alcohols
UOM	1 EACH

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79289>