

Semiprep/Preparative Guard Column, Carbonate Form

Specification Sheet (Part/REF # 79919)



HC-75 is a unique cation exchange resin designed specifically for separating and purifying oligosaccharides with a degree of polymerization ≤ 8 . Its specialized functionality makes it an essential tool in food science, biochemistry, and pharmaceutical manufacturing, where detailed carbohydrate analysis is critical.

Product Specifications

Separation Mode	Cation Exchange
Stationary Phase	PRP-X500
Format	Guard Column
Hardware Inner Diameter	4.6 mm
Hardware Length	33 mm
Particle Size	12-20 μm
Packing Material Type	PSDVB/Sulfonic Acid
Particle Surface Area (m ² /cm ³)	1.4
Pore Volume (cm ³ /g)	0.0063
Guard Column Size	Semiprep/Preparative
Hardware Material	Stainless Steel
Max Pressure	400 psi
UOM	1 EACH

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79919>