

PRP-3 Analytical Guard Column

Specification Sheet (Part/REF # 79295)



Organic compounds: large molecules (> 2,000 mw), peptides, proteins, protein digests, protected and deprotected oligonucleotides, nucleic acids.

Product Specifications

Separation Mode	Reversed-Phase
Stationary Phase	PRP-3
Format	Guard Column
Hardware Inner Diameter	2.1 mm
Hardware Length	33 mm
Particle Size	10 µm
Packing Material Type	PSDVB
Particle Surface Area (m ² /cm ³)	48
Pore Volume (cm ³ /g)	0.20
Guard Column Size	Analytical
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Compound Classification	Organics
UOM	1 EACH

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79295>