

PRP-X110 Anion Exchange Analytical Guard Column

Specification Sheet (Part/REF # 79293)



High sensitivity detection ion exchange chromatography of organic and inorganic anions, organic acids, organic and inorganic arsenic species, halides, nucleotides.

Product Specifications

Separation Mode	Anion Exchange
Stationary Phase	PRP-X110
Format	Guard Column
Hardware Inner Diameter	2.1 mm
Hardware Length	33 mm
Particle Size	12-20 µm
Packing Material Type	PSDVB/Trimethylammonium
Particle Surface Area (m2/cm3)	476
Pore Volume (cm3/g)	0.728
Guard Column Size	Analytical
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Compound Classification	Organic/Inorganic Anions
UOM	1 EACH

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/hplc-columns/79293>