

0.5 mm, TLL-Male Venting Needle, 51 mm, point style 5, 1/PK

Specification Sheet (Part/REF # 210544)



Hamilton venting needles are made of high grade stainless steel and are ideally suited for a wide range of applications. Various lengths and different hub connections are available to perfectly match the requirements of laboratory equipment manufacturers.

Product Specifications

Gauge	0.25-0.50 mm
Point Style	5: Conical Side Port
Needle Length	2 inch (51 mm)
Autoclavable	No
Needle Inner Diameter	0.010 inches (0.25 mm)
Needle Outer Diameter	0.020 inches (0.50 mm)
Sales Quantity	1 needle per pack
Hub	TLL-Male

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/needles/210540-1>