

# 1.0 mm, TLL-Male Venting Needle, 76 mm, point style 5, 1/PK

Specification Sheet (Part/REF # 210545)



Hamilton venting needles are made of high grade stainless steel and are ideally suited for a wide range of applications. Various lengths and different hub connections are available to perfectly match the requirements of laboratory equipment manufacturers.

## Product Specifications

Gauge	0.70-1.00 mm
Point Style	5: Conical Side Port
Needle Length	3 inch (76 mm)
Autoclavable	No
Needle Inner Diameter	0.028 inches (0.70 mm)
Needle Outer Diameter	0.039 inches (1.00 mm)
Sales Quantity	1 needle per pack
Hub	TLL-Male

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/needles/210545>