

# 1.6 mm, C-Male Venting Needle, 102 mm, point style 5, 1/PK

Specification Sheet (Part/REF # 210554)



Hamilton venting needles are made of high grade stainless steel and are ideally suited for a wide range of applications. Various lengths and different hub connections are available to perfectly match the requirements of laboratory equipment manufacturers.

## Product Specifications

Gauge	1.20-1.60 mm
Point Style	5: Conical Side Port
Needle Length	4.00 in (101.6 mm)
Autoclavable	Yes
Needle Inner Diameter	0.047 inches (1.20 mm)
Needle Outer Diameter	0.063 inches (1.60 mm)
Sales Quantity	1 needle per pack
Hub	C-Male

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/needles/210554>