

2.0 mm, TLL-Female Venting Needle, 102 mm, point style 5, 1/PK

Specification Sheet (Part/REF # 210543)



Hamilton venting needles are made of high grade stainless steel and are ideally suited for a wide range of applications. Various lengths and different hub connections are available to perfectly match the requirements of laboratory equipment manufacturers.

Product Specifications

Point Style	5: Conical Side Port
Needle Length	4.00 in (101.6 mm)
Autoclavable	Yes
Needle Inner Diameter	0.063 inches (1.60 mm)
Needle Outer Diameter	0.079 inches (2.00 mm)
Sales Quantity	1 needle per pack
Hub	TLL-Female

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/needles/210543>