

Priming Kit

Specification Sheet (Part/REF # PRMKIT)



The Priming Kit is used to backfill syringes to eliminate compressible air from the system. This is especially useful when the syringe is connected to tubing that exceeds the volume of the syringe.

Product Specifications

Consists of	one 250 uL syringe (p/n 81120), one 6 pack of 30 gauge needles, and one pack of septa (p/n 75826).
Compatible Syringes	The priming kit is compatible with 600, 700, 1700, removable needle syringes and 7000 series syringes

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/syringes/prmkita>

Spec. Version B