

10 µL Conductive Sterile Tips

Specification Sheet (Part/REF # 235935)



10 µL Conductive Sterile Non-Filter Tips

Product Specifications

Platform Compatibility	Microlab VANTAGE Microlab STAR Microlab NIMBUS
Tip Type	Standard
Tip Volume	10 µL
Maximum Tip Capacity	10 µL
Filter	Non-Filter
Sterile	Sterile
Tip Color (Conductivity)	Black (Conductive)
Barcode Label Color	Orange
Tip Material	Polypropylene - PP
Tip Recycling Code	5
Rack Material	Polypropylene - PP
Rack Recycling Code	5
Pipetting Compatibility	1000 µL channel - 96 MultiProbe Head
Purity Level	Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, and PCR inhibitors
Tips / Rack	96 Tips/Rack
SLAS ANSI	N/A
Tip Length	29.9 mm
Orifice Size	0.3 mm
Case Size	Case of 5,760 tips (Blister 5 x 96 tips per rack)
Packaging Option	Blister 96 Tip

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/disposable-tips/10-µl-conductive-sterile-non-filter-tips>

Spec. Version B