


300 µL Conductive Filter Tips

Specification Sheet (Part/REF # 235903)



300 µL Conductive Non-Sterile Filter Tips

Product Specifications

Platform Compatibility	Microlab VANTAGE Microlab STAR Microlab NIMBUS Microlab Prep
Tip Type	Standard
Tip Volume	300 µL
Maximum Tip Capacity	300 µL
Filter	Filter
Sterile	Non-Sterile
Tip Color (Conductivity)	Black (Conductive)
Barcode Label Color	Yellow
Tip Material	Polypropylene - PP
Tip Recycling Code	5
Rack Material	Polypropylene - PP
Rack Recycling Code	5
Pipetting Compatibility	1000 µL channel - 96 MultiProbe Head
Purity Level	Free of DNase, RNase, Pyrogen (endotoxins), Human DNA, and PCR inhibitors
Tips / Rack	96 Tips/Rack
SLAS ANSI	N/A
Tip Length	59.9 mm
Orifice Size	0.4 mm
Case Size	Case of 5,760 tips (Blister 5 x 96 tips per rack)
Packaging Option	Blister 96 Tip

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/disposable-tips/300-µl-conductive-non-sterile-filter-tips>

Spec. Version B