

# ClickSure Tips 1000 µL, Bulk, Standard, 960 tips/bag

Specification Sheet (Part/REF # 235545)



The Disposable Tip Hand Probe utilizes Hamilton's unique ClickSure® tips which have been validated for superior liquid handling.

## Product Specifications

|                       |               |
|-----------------------|---------------|
| Disposable Tip Volume | 1000 µL       |
| Packaging             | Bulk          |
| Material              | Polypropylene |

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/pipettes/235545>