

50 µL Nested Conductive Tips

Specification Sheet (Part/REF # 235947)



50 µL Nested Conductive Non-Sterile Non-Filter Tips

Product Specifications

Platform Compatibility	Microlab VANTAGE Microlab STAR Microlab NIMBUS Microlab Prep
Tip Type	Standard
Tip Volume	50 µL
Maximum Tip Capacity	50 µL
Filter	Non-Filter
Sterile	Non-Sterile
Tip Color (Conductivity)	Black (Conductive)
Barcode Label Color	Light Red
Tip Material	Polypropylene - PP
Tip Recycling Code	5
Rack Material	Polystyrene - PS
Rack Recycling Code	6
Pipetting Compatibility	1000 µL channel - 96 MultiProbe Head
Tips / Rack	96 Tips/Rack
Tip Length	50.4 mm
Orifice Size	0.4 mm
Case Size	Case of 11,520 tips (NTR 5 x 4 stack)
Packaging Option	Nested 96 Tip

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/disposable-tips/50-ul-nested-conductive-tips-1-1>

Spec. Version B