


# Septa, Long Life GC Autosampler, 11.0 mm DIA (25/PK)

Specification Sheet (Part/REF # 75826)



Long Life GC Autosampler Septa are designed for use with GC Autosamplers. The 5, 9, 11 and 11.5 mm diameter Long Life septa are available with a unique pre-pierced center depression to guide the syringe needle to the same injection point, extending septa life, reducing needle bending and eliminates the common type of coring that occurs on the first injection.

## Product Specifications

Diameter	11 mm (0.433 inch)
Maximum Temperature	400°C / 752°F
Thickness	3.00 mm (0.118 inch)
Sales Quantity	25/pack

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/gc-septa/75826>

Spec. Version B