

# Semiprep/Preparative Guard Cartridge Holder, Stainless Steel

Specification Sheet (Part/REF # 5095-01)



PRP-X400 is a type of cation exchange resin, highly effective in detecting and removing the pesticide glyphosate and its metabolites from drinking water. In addition to this specific application, it handles a broad range of inorganic and organic cations, making it indispensable in numerous scientific and industrial contexts.

## Product Specifications

Format	Guard Cartridge Holder
Guard Column Size	Semiprep/Preparative
Hardware Material	Stainless Steel
UOM	1 EACH

Specifications are subject to change without notice  
<https://www.hamiltoncompany.com/hplc-columns/5095-01>

Spec. Version B