

PRP-X100 Anion Exchange Semiprep/Preparative Guard Cartridge Replacement Cartridges (2/pk), Stainless Steel

Specification Sheet (Part/REF # 79126)



High capacity ion exchange chromatography of organic and inorganic anions, organic acids, organic and inorganic arsenic species, halides, nucleotides.

Product Specifications

Separation Mode	Anion Exchange
Stationary Phase	PRP-X100
Format	Replacement Guard Cartridge
Hardware Inner Diameter	4.6 mm
Hardware Length	20 mm
Particle Size	12-20 µm
Packing Material Type	PSDVB/Trimethylammonium
Particle Surface Area (m ² /cm ³)	405
Pore Volume (cm ³ /g)	0.64
Guard Column Size	Semiprep/Preparative
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Organic/Inorganic Anions
UOM	2/pk

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79126>

Spec. Version B