

RCX-10 Analytical Guard Column Starter Kit (1 holder, 2 cartridges), Stainless Steel

Specification Sheet (Part/REF # 79462)



High exchange capacity anion exchange column for carbohydrate separation: Carbohydrates, polysaccharides, sugar oligomers up to DP8.

Product Specifications

Separation Mode	Anion Exchange
Stationary Phase	RCX-10
Format	Guard Cartridge Kit
Hardware Inner Diameter	2.1 mm
Hardware Length	20 mm
Particle Size	12-20 μm
Packing Material Type	PSDVB/Trimethylammonium
Particle Surface Area (m ² /cm ³)	358
Pore Volume (cm ³ /g)	0.624
Guard Column Size	Semiprep/Preparative
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Carbohydrates
UOM	1 holder, 2 cartridges

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79462>

Spec. Version B