

RCX-10 Analytical Guard Column Replacement Cartridges (5/pk), Stainless Steel

Specification Sheet (Part/REF # 79463)



High exchange capacity anion exchange column for carbohydrate separation: Carbohydrates, polysaccharides, sugar oligomers up to DP8.

Product Specifications

Separation Mode	Anion Exchange
Stationary Phase	RCX-10
Format	Replacement Guard Cartridge
Hardware Inner Diameter	2.1 mm
Hardware Length	20 mm
Particle Size	12-20 µm
Packing Material Type	PSDVB/Trimethylammonium
Particle Surface Area (m2/cm3)	358
Pore Volume (cm3/g)	0.624
Guard Column Size	Analytical
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Carbohydrates
UOM	5/pk

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79463>

Spec. Version B