

RCX-30, 12-20 µm Bulk Resin

Specification Sheet (Part/REF # 79706)



Increased exchange capacity anion exchange columns for better separation of carbohydrate mixtures: Complex mono- and disaccharide separation.

Product Specifications

Separation Mode	Anion Exchange
Stationary Phase	RCX-30
Format	Bulk Resin
Particle Size	12-20 µm
Packing Material Type	PSDVB/Trimethylammonium
Pore Size	100 Å
Particle Surface Area (m²/cm³)	105
Pore Volume (cm³/g)	0.28
Exchange Capacity	1.0 meq/gm
Restoration	Flush with 150 mL 0.1 N sodium hydroxide.
UOM	per gram

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79706>

Spec. Version B