

PRP-X800 Analytical Guard Column Starter Kit (1 holder, 2 cartridges), PEEK

Specification Sheet (Part/REF # 79831)



Protect your PRP-X800 column.

Product Specifications

Separation Mode	Cation Exchange
Stationary Phase	PRP-X800
Format	Guard Cartridge Kit
Hardware Inner Diameter	3.0 mm
Hardware Length	8 mm
Particle Size	7 μm
Packing Material Type	PSDVB/Itaconate
Particle Surface Area (m ² /cm ³)	50
Pore Volume (cm ³ /g)	0.166
Guard Column Size	Analytical
Hardware Material	PEEK
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Mono-Divalent Organic/Inorganic Cations - Transition Metals
Frit Material	Titanium
UOM	1 holder, 2 cartridges

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79831>

Spec. Version B