

PRP-X300 Analytical Guard Cartridge Replacement Cartridges (5/pk), PEEK

Specification Sheet (Part/REF # 79374)



Protect your PRP-X300 column.

Product Specifications

Separation Mode	Ion Exclusion
Stationary Phase	PRP-X300
Format	Replacement Guard Cartridge
Hardware Inner Diameter	3.0 mm
Hardware Length	8 mm
Particle Size	12-20 μm
Packing Material Type	PSDVB/Sulfonic Acid
Particle Surface Area (m ² /cm ³)	360
Pore Volume (cm ³ /g)	0.748
Guard Column Size	Analytical
Hardware Material	PEEK
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Organic Acids/Alcohols
Frit Material	Titanium
UOM	5/pk

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79374>

Spec. Version B