

PRP-X300 Analytical Guard Cartridge Starter Kit (1 holder, 2 cartridges), Stainless Steel

Specification Sheet (Part/REF # 79460)



Protect your PRP-X300 column.

Product Specifications

Separation Mode	Ion Exclusion
Stationary Phase	PRP-X300
Format	Guard Cartridge Kit
Hardware Inner Diameter	2.1 mm
Hardware Length	20 mm
Particle Size	12-20 μm
Packing Material Type	PSDVB/Sulfonic Acid
Particle Surface Area (m ² /cm ³)	360
Pore Volume (cm ³ /g)	0.748
Guard Column Size	Semiprep/Preparative
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Organic Acids/Alcohols
UOM	1 holder, 2 cartridges

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79460>

Spec. Version B