

PRP-X300 Semiprep/Preparative Guard Cartridge Replacement Cartridges (2/pk), Stainless Steel

Specification Sheet (Part/REF # 79130)



Protect your PRP-X300 column.

Product Specifications

Separation Mode	Ion Exclusion
Stationary Phase	PRP-X300
Format	Replacement Guard Cartridge
Hardware Inner Diameter	4.6 mm
Hardware Length	20 mm
Particle Size	12-20 µm
Packing Material Type	PSDVB/Sulfonic Acid
Particle Surface Area (m ² /cm ³)	360
Pore Volume (cm ³ /g)	0.748
Guard Column Size	Semiprep/Preparative
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Exchange Capacity	0.17 meq/gm
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Temp Limits	5-60°C
Compound Classification	Organic Acids/Alcohols
UOM	2/pk

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79130>

Spec. Version B