

PRP-1 Semiprep/Preparative Guard Column Starter Kit (1 holder, 1 cartridge), Stainless Steel

Specification Sheet (Part/REF # 79121)



Organic compounds: small molecules (< 2,000 mw), oligonucleotides, pharmaceuticals, steroids, nucleic acids, vitamins, herbicides.

Product Specifications

Separation Mode	Reversed-Phase
Stationary Phase	PRP-1
Format	Guard Cartridge Kit
Hardware Inner Diameter	4.6 mm
Hardware Length	20 mm
Particle Size	12-20 µm
Packing Material Type	PSDVB
Particle Surface Area (m ² /cm ³)	400
Pore Volume (cm ³ /g)	0.772
Guard Column Size	Semiprep/Preparative
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Organics
UOM	1 holder, 1 cartridge

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79121>

Spec. Version B