

PRP-1 Analytical Guard Cartridge Starter Kit (1 holder, 2 cartridges), Stainless Steel

Specification Sheet (Part/REF # 79447)



Organic compounds: small molecules (< 2,000 mw), oligonucleotides, pharmaceuticals, steroids, nucleic acids, vitamins, herbicides.

Product Specifications

Separation Mode	Reversed-Phase
Stationary Phase	PRP-1
Format	Guard Cartridge Kit
Hardware Inner Diameter	2.1 mm
Hardware Length	20 mm
Particle Size	10 μ m
Packing Material Type	PSDVB
Particle Surface Area (m ² /cm ³)	400
Pore Volume (cm ³ /g)	0.772
Guard Column Size	Analytical
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Organics
UOM	1 holder, 2 cartridges

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79447>

Spec. Version B