

PRP-3 Analytical Guard Column Starter Kit (1 holder, 2 cartridges), PEEK

Specification Sheet (Part/REF # 79393)



Protect your PRP-3 column.

Product Specifications

Separation Mode	Reversed-Phase
Stationary Phase	PRP-3
Format	Guard Cartridge Kit
Hardware Inner Diameter	3.0 mm
Hardware Length	8 mm
Particle Size	10 μ m
Packing Material Type	PSDVB
Particle Surface Area (m ² /cm ³)	48
Pore Volume (cm ³ /g)	0.20
Guard Column Size	Analytical
Hardware Material	PEEK
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Peptides/Proteins
Frit Material	Titanium
UOM	1 holder, 2 cartridges

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79393>

Spec. Version B