

PRP-3 Semiprep/Preparative Guard Column Replacement Cartridges (2/pk), Stainless Steel

Specification Sheet (Part/REF # 79124)



Protect your PRP-3 column.

Product Specifications

Separation Mode	Reversed-Phase
Stationary Phase	PRP-3
Format	Replacement Guard Cartridge
Hardware Inner Diameter	4.6 mm
Hardware Length	20 mm
Particle Size	12-20 μm
Packing Material Type	PSDVB
Particle Surface Area (m ² /cm ³)	48
Pore Volume (cm ³ /g)	0.20
Guard Column Size	Semiprep/Preparative
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Peptides/Proteins
UOM	2/pk

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79124>

Spec. Version B