

PRP-C18 Analytical Guard Column Replacement Cartridges (5/pk), Stainless Steel

Specification Sheet (Part/REF # 79686)



Protect your PRP-C18 column.

Product Specifications

Separation Mode	Reversed-Phase
Stationary Phase	PRP-C18
Format	Replacement Guard Cartridge
Hardware Inner Diameter	2.1 mm
Hardware Length	20 mm
Particle Size	12-20 µm
Packing Material Type	PS-DVB functionalized with C18
Particle Surface Area (m2/cm3)	368
Pore Volume (cm3/g)	0.624
Guard Column Size	Analytical
Hardware Material	Stainless Steel
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Organics
UOM	5/pk

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/hplc-columns/79686>

Spec. Version B