

PRP-C18 Analytical Guard Column Starter Kit (1 holder, 2 cartridges), PEEK

Specification Sheet (Part/REF # 79687)



Protect your PRP-C18 column.

Product Specifications

Separation Mode	Reversed-Phase
Stationary Phase	PRP-C18
Format	Guard Cartridge Kit
Hardware Inner Diameter	3.0 mm
Hardware Length	8 mm
Particle Size	12-20 μm
Packing Material Type	PS-DVB functionalized with C18
Particle Surface Area (m^2/cm^3)	368
Pore Volume (cm^3/g)	0.624
Guard Column Size	Analytical
Hardware Material	PEEK
Max Pressure	5,000 psi
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Compound Classification	Organics
Frit Material	Titanium
UOM	1 holder, 2 cartridges

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79687>

Spec. Version B