

# PRP-X100 Anion Exchange HPLC Column, PEEK, 2.1 x 250 mm, 5 µm

Specification Sheet (Part/REF # 79670)



Organic and inorganic anions, and arsenic speciation.

## Product Specifications

Separation Mode	Anion Exchange
Stationary Phase	PRP-X100
Format	HPLC Column
Hardware Inner Diameter	2.1 mm
Hardware Length	250 mm
Particle Size	5 µm
Packing Material Type	PSDVB/Trimethylammonium
Particle Surface Area (m <sup>2</sup> /cm <sup>3</sup> )	405
Pore Volume (cm <sup>3</sup> /g)	0.64
Hardware Material	PEEK
Max Pressure	5,000 psi
Exchange Capacity	0.19 meq/gm
Exchange Capacity/Column	0.076 meq/column
Mobile Phase Limits	pH 1-13. 0-100% aqueous, organic modifier
Temp Limits	5-30°C
Restoration	Pump approximately 50 mL of methanol with 1% 6 N nitric acid.
Compound Classification	Organic/Inorganic Anions

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79670>

Spec. Version B