

HC-75 (Ca²⁺), 9 µm Bulk Resin

Specification Sheet (Part/REF # 79709)



Pack your own column.

Product Specifications

Separation Mode	Cation Exchange
Stationary Phase	HC-75 (Ca ²⁺)
Format	Bulk Resin
Particle Size	9 µm
USP Classification	USP L19
Pore Size	Gel-Type
Exchange Capacity	5 meq/gm
Restoration	Flush with 1% calcium chloride at 0.1 mL/min overnight.
UOM	per gram

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79709>

Spec. Version B