

HC-75 (Pb 2+), 9 µm Bulk Resin

Specification Sheet (Part/REF # 79712)



Pack your own column.

Product Specifications

Separation Mode	Cation Exchange
Stationary Phase	HC-75 (Pb+2)
Format	Bulk Resin
Particle Size	9 µm
USP Classification	USP L34
Pore Size	Gel-Type
Exchange Capacity	5 meq/gm
Restoration	Flush with 1% lead nitrate at 0.1 mL/min overnight.
UOM	per gram

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/hplc-columns/79712>

Spec. Version B