

PRP-3, 10 µm Bulk Resin

Specification Sheet (Part/REF # 79701)



Pack your own column.

Product Specifications

Separation Mode	Reversed-Phase
Stationary Phase	PRP-3
Format	Bulk Resin
Particle Size	10 µm
USP Classification	USP L21
Packing Material Type	PSDVB
Pore Size	300 Å
Particle Surface Area (m²/cm³)	48
Pore Volume (cm³/g)	0.20
Max Pressure	5,000 psi
Restoration	Run a gradient of 100% water to 100% acetonitrile (Repeat three times)
UOM	per gram

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/hplc-columns/79701>

Spec. Version B